

AmpliSeq™ for Illumina BRCA Panel

Fast and accurate detection of somatic and germline mutations in *BRCA1* and *BRCA2*.

Highlights

- **Relevant gene content**
Target all exonic regions and flanking intronic sequences of *BRCA1* and *BRCA2*
- **Fast, streamlined workflow**
Prepare sequencing-ready libraries in a single day using as little as 1 ng high-quality DNA or 10 ng DNA from FFPE tissue
- **Accurate data**
Detect germline and somatic mutations down to 5% frequency using local or cloud-based analysis

Introduction

The AmpliSeq for Illumina *BRCA* Panel is a targeted resequencing assay designed for detecting somatic and germline mutations across all exonic regions and the flanking intronic sequences of *BRCA1* and *BRCA2* (Table 1). *BRCA1* and *BRCA2* are human tumor suppressor genes that, when carrying specific mutations, have been implicated in an increased risk for breast and ovarian cancers.¹ Understanding BRCA status within the tumor may be a factor when researching potential therapies.²

To assist with quick and accurate assessment of genomic variation within *BRCA1* and *BRCA2*, Illumina offers the AmpliSeq for Illumina *BRCA* Panel. The *BRCA* panel is part of a streamlined workflow that includes PCR-based library preparation, Illumina sequencing by synthesis (SBS) chemistry and next-generation sequencing (NGS) technology, and automated analysis. Requiring as little as 1 ng high-quality DNA per pool, the two-pool panel can be used with low-quality and low-quantity samples, including formalin-fixed, paraffin-embedded (FFPE) tissues.

Simple, streamlined workflow

The AmpliSeq for Illumina *BRCA* Panel is part of a DNA-to-variant solution that offers streamlined content, easy-to-perform library preparation, push-button sequencing systems, and simplified data analysis.

Library preparation follows a straightforward, PCR-based protocol that can be completed in as little as 5 hours, with < 1.5 hours hands-on time. Resulting libraries can be normalized, pooled, and then loaded on to a flow cell for sequencing. Prepared libraries are sequenced using proven SBS chemistry on any compatible Illumina sequencing system (Table 2).

Resulting data can be analyzed locally with Local Run Manager or easily streamed into BaseSpace™ Sequence Hub. Local Run Manager and BaseSpace Sequence Hub can access the DNA Amplicon analysis workflow to perform alignment and variant

calling. BaseSpace Sequence Hub provides access to BaseSpace Variant Interpreter, which assists with turning variant call data into annotated results.

Table 1: AmpliSeq for Illumina BRCA Panel at a glance

| Parameter | Specification |
|--|--|
| No. of genes | 2 |
| Targets | All exonic regions of the <i>BRCA1</i> and <i>BRCA2</i> tumor suppressor genes and flanking intronic sequences |
| Cumulative target size | 22 kb |
| Variant types | SNVs, indels ^a |
| Amplicon size | 98 bp on average |
| No. of amplicons | 265 |
| Input DNA requirement | 1-100 ng (10 ng recommended per pool) |
| No. of pools per panel | 2 |
| Supported sample types | FFPE tissue, blood |
| <i>Germline</i> : percent targets covered at minimum 50× at recommended throughput | > 95% |
| <i>Somatic</i> : percent targets covered at minimum 500× at recommended throughput | > 95% |
| Coverage uniformity (percent of targets with >0.2× mean coverage) | > 95% |
| Percent on-target aligned reads | > 90% |
| Total assay time ^b | 5 hours |
| Hands-on time | < 1.5 hours |
| DNA-to-data time | 2.5 days |

a. SNVs: single nucleotide variations; indels: insertions/deletions
b. Time represents library preparation only and does not include library quantification, normalization, or pooling.

Data on file at Illumina, Inc. 2017

Table 2: Illumina sequencing systems recommended for use with the AmpliSeq for Illumina BRCA Panel

| Instrument | No. of Samples per Run | Run Time |
|---|------------------------|----------|
| For investigating somatic mutations | | |
| iSeq™ 100 System | 12 | 17 hours |
| MiniSeq™ System (mid output) | 24 | 17 hours |
| MiniSeq System (high output) | 80 | 24 hours |
| MiSeq System (v2 chemistry Nano) | 3 | 17 hours |
| MiSeq System (v2 chemistry Micro) | 12 | 19 hours |
| MiSeq System (v2 chemistry) | 48 | 24 hours |
| MiSeq System (v3 chemistry) | 80 | 32 hours |
| For investigating germline mutations | | |
| iSeq 100 System | 96 | 17 hours |
| MiSeq System (v2 chemistry Nano) | 32 | 17 hours |
| MiSeq System (v2 chemistry Micro) | 96 | 19 hours |

A maximum of 96 samples can be multiplexed.




Figure 1: High Coverage Uniformity and On-Target Alignment—DNA extracted from HD and breast tissue samples of varying quality was evaluated using the AmpliSeq for Illumina *BRCA* Panel and sequenced on the MiSeq System. Error bars indicate variability of technical replicates. ΔCq is an indicator of the quality of the DNA.

 [Learn more about AmpliSeq for Illumina Informatics Solutions >](#)

Accurate data

To demonstrate assay capabilities, Horizon Discovery (HD) and breast tissue (tumor/normal pair) samples were evaluated using the AmpliSeq for Illumina *BRCA* Panel on the MiSeq™ System. Results showed high coverage uniformity and on-target percentage of aligned reads, even with varying sample quality (Figure 1). In addition, two HD samples of varying quality were evaluated for variant calling accuracy. Data showed high concordance between expected and detected SNVs (Table 3).

Ordering information

Order AmpliSeq for Illumina products online at www.illumina.com

| Product | Catalog No. |
|---|-------------|
| AmpliSeq for Illumina <i>BRCA</i> Panel (24 reactions) | 20019168 |
| AmpliSeq for Illumina Library PLUS (24 reactions) | 20019101 |
| AmpliSeq for Illumina Library PLUS (96 reactions) | 20019102 |
| AmpliSeq for Illumina Library PLUS (384 reactions) | 20019103 |
| AmpliSeq for Illumina CD Indexes Set A (96 indexes, 96 samples) | 20019105 |
| AmpliSeq for Illumina Direct FFPE DNA | 20023378 |
| AmpliSeq for Illumina Library Equalizer | 20019171 |

Learn more

Learn more about the [AmpliSeq for Illumina *BRCA* Panel](#)

Learn more about the [AmpliSeq for Illumina targeted sequencing solution](#)

References

- BRCA1 and BRCA2: Cancer Risk and Genetic Testing Fact Sheet - National Cancer Institute. <https://www.cancer.gov/about-cancer/causes-prevention/genetics/brca-fact-sheet#q1>. Accessed October 30, 2017.

- Approved Drugs > FDA approves olaparib tablets for maintenance treatment in ovarian cancer. <https://www.fda.gov/drugs/informationondrugs/approved-drugs/ucm572143.htm>. Accessed November 21, 2017.

Table 3: High Concordance Between Expected and Detected Variant Frequency

| Known Variant | Expected VF ^a | Detected VF ^a | No. of Samples | Call Rate |
|--|--------------------------|--------------------------|----------------|-----------|
| HD795 (Horizon Discovery high-quality sample) | | | | |
| <i>BRCA2</i> N289H | 7.5% | 7.4% | 16 | 100% |
| <i>BRCA2</i> N991D | 7.5% | 7.2% | 16 | 100% |
| <i>BRCA2</i> D1420Y | 32.5% | 32.5% | 16 | 100% |
| <i>BRCA2</i> V2466A | 100% | 99.9% | 16 | 100% |
| <i>BRCA1</i> S1613G | 7.5% | 7.4% | 16 | 100% |
| <i>BRCA1</i> R1443STOP | 32.5% | 32.8% | 16 | 100% |
| <i>BRCA1</i> K1183R | 7.5% | 7.9% | 16 | 100% |
| <i>BRCA1</i> P871L | 15.0% | 14.5% | 16 | 100% |
| <i>BRCA1</i> K820E | 7.5% | 7.7% | 16 | 100% |
| <i>BRCA1</i> D435Y | 7.5% | 6.9% | 16 | 100% |
| HD810 (Horizon Discovery formalin-fixed sample) | | | | |
| <i>BRCA2</i> N289H | 7.5% | 7.0% | 10 | 100% |
| <i>BRCA2</i> N991D | 7.5% | 7.6% | 10 | 100% |
| <i>BRCA2</i> D1420Y | 32.5% | 33.9% | 10 | 100% |
| <i>BRCA2</i> V2466A | 100% | 99.9% | 10 | 100% |
| <i>BRCA1</i> S1613G | 7.5% | 7.3% | 10 | 100% |
| <i>BRCA1</i> R1443STOP | 32.5% | 33.9% | 10 | 100% |
| <i>BRCA1</i> K1183R | 7.5% | 7.5% | 10 | 100% |
| <i>BRCA1</i> P871L | 15.0% | 14.6% | 10 | 100% |
| <i>BRCA1</i> K820E | 7.5% | 7.7% | 10 | 100% |
| <i>BRCA1</i> D435Y | 7.5% | 7.1% | 10 | 100% |

a. VF: variant frequency

DNA from high-quality (HD795) and formalin-fixed (HD810) samples was evaluated using the AmpliSeq for Illumina *BRCA* Panel and sequenced on the MiSeq System. Results show that 100% of expected SNVs were detected.